

Uroflowmeters and uroflowmetric transducers

DanFlow 1000 and 1100

Compact systems for measurement control and evaluation of signal from flow sensors – uroflowmeters.



- They include an inbuilt thermal printer and unique automatic archiving of measurement data onto an SD card.

- DanFlow 1000 connects uroflowmetric transducer by cable.

- DanFlow 1100 connects uroflowmetric transducer wirelessly.



Weight uroflowmetric transducers, two height adjustable stands, commode, transducer can be fixed to diagnostic couch TryTable.



DanFlow 2000 connects by a cable to a standard computer using a Windows O/S.



DanFlow 3000 connects wirelessly to a standard computer fitted with a patented USB radio adapter.

Accessories



External EMG amplifier



Uroflowmeter with examination chair



Uroflowmeter on economical stand



Profilometric puller on a stand

