





Versatile endovaginal transducer

| 8819 |



- Superior image quality for all-round gynecology examinations
- Simple and convenient needle guide
- Ergonomic and slim design



Endovaginal cross-section of uterus showing enlarged endometrium.



Endovaginal scan of uterus showing large suspicious mass.

For your routine gynecology investigations, interventional gynecology and IVF*.

■ Superior image quality for all-round gynecology examinations

The 8819 endovaginal endfire transducer is the ideal all-round ultrasound tool for routine gynecological investigations and IVF.

The 8819 offers good penetration and gives you detailed images at all depths.

■ Simple and convenient needle guide

The 8819 needle guide clicks easily and securely into place and helps simplify interventional procedures.

■ Ergonomic and slim design

Ergonomic and slim design for comfortable use. The 8819's long, slim design minimizes the discomfort of a transvaginal examination.



The 8819 needle guide helps simplify intervention.



■ World Headquarters:

Mileparken 34 · DK-2730 · Herlev · Denmark

Tel: (45) 4452 8100 · Fax: (45) 4452 8199 · www.bkmed.com

A wholly owned subsidiary of Analogic Corporation

Specifications

Frequency range

5-9 MHz

Penetration

11 cm

Sector angle

150°

Contact surface

26 x 5 mm

Physical data

22 x 35 x 319 mm

Weight

140 g

**8819 is not cleared for IVF use in the USA.*