




## The New Gold Standard in Neurosurgical Ultrasound



*Easy-to-use sterile needle guides*



*World's smallest high-resolution burr hole transducer*



*High-resolution craniotomy transducer*

### Real Clinical Impact with High-Resolution Neurosurgical Imaging

- Premium Performance Neurosurgical Imaging
- Exceptionally Intuitive System
- Unmatched Transducer Function and Design
- Convenient, Sterile Needle Guides



# Real Clinical Impact with High-Resolution Neurosurgical Imaging



Large brain tumor evaluated during surgery using craniotomy transducer



Parasitic brain cyst in coronal plane – using craniotomy transducer



Parasitic brain cyst in sagittal plane – using craniotomy transducer



Intraoperative imaging of spinal tumor in transverse plane using craniotomy transducer

## Premium Performance Neurosurgical Imaging

- Superior, detailed images for precise neurostructural imaging
- Clear, sharp images of brain and spinal cord
- Real-time image guidance

## Exceptionally Intuitive System

- Clinically intuitive – easy-to-use
- Small system footprint for added OR convenience
- Up to 4 hours of plug-free imaging

## Unmatched Transducer Function and Design

- World's smallest high-resolution burr-hole transducer
- High-resolution craniotomy transducer
- Compatible with modern sterilization techniques
- Direct contact for sharper images

## Convenient, Sterile Needle Guides

- Easy to use – needle guide just snaps on and off
- One patient, one guide – sterile and ready-to-use
- Easily detachable – needle and shunt stay in place

flexFocus 800



Quantum Technology



Specifications	Frequency Range	Focal Range (typical)	Contact Surface (overall)	Sector Angle
Craniotomy Transducer – 8862	8862: 10–3.8 MHz	8862: 5–68 mm	8862: 29 x10 mm	8862: 66°
Burr Hole Transducer – 8863	8863: 10–3.8 MHz	8863: 5–57 mm	8863: 10 x 8.6 mm	8863: 46°

**Analogic Corporation - Headquarters USA**  
 8 Centennial Drive, Peabody, MA 01960  
 T 978-326-4000  
 analogic.com

**BK Medical - Sales and Service USA**  
 8 Centennial Drive, Peabody, MA 01960  
 T 978-326-1300  
 bkmed.com

**BK Medical - Europe and Rest of World**  
 Mileparken 34, 2730, Herlev, Denmark  
 T +45 4452 8100  
 F +45 4452 8199  
 bkmed.com



**BK Medical ApS, Mileparken 34, 2730 Herlev, Denmark. Tel.: +45 44528100 Fax: +45 44528199 Email: info@bkmed.dk**  
 The name Analogic is a registered trademark of Analogic Corporation. The globe logo is a trademark of Analogic Corporation.